



WellDone!

Plate setup and PCR setup for the Roche qPCR systems

welldone-prompter.github.io/WellDone

Free & offline

Windows · one single file · no installation
no administrator rights · German & English

WORKFLOW From sample entry to the finished plate



What WellDone! does

- ▶ **Record samples** – manually, with a barcode scanner or as an import from MagNA Pure 96/24 and Excel
- ▶ **Fill the plate automatically** – 96 and 384 wells, including standards and controls
- ▶ **Assay library** – set up reagents, controls and PCR profiles once and reuse them
- ▶ **MasterMix calculation** – required volumes including a safety margin
- ▶ **Pipetting scheme as a PDF** – ready to print and colour-matched to the plate
- ▶ **Export** – setup files for LightCycler PRO, 480 and 96 plus the order for MagNA Pure 96
- ▶ **Bilingual** – interface and manual in German and English

WellDone! — Pipettierschema | 260627_Gastro_01
Erstellt: 27.06.2025 22:21

1. MasterMix ansetzen und je Well verteilen

Komponente (µl)	Gastro Virus je 1. Reil	Gastro Virus present (25x)	Gastro Bacteria je 1. Reil	Gastro Bacteria present (25x)	Gastro Ecoli je 1. Reil	Gastro Ecoli present (25x)
Wasser	9.4	235.0	9.9	247.5	9.4	235.0
Primer/Probes	1.5	37.5	1.9	47.5	1.5	37.5
MasterMix	4.0	100.0	5.0	200.0	4.0	100.0
RT-Enzym	0.1	2.5	0.2	5.0	0.1	2.5
Summe MasterMix	15.0	375.0	20.0	500.0	15.0	375.0
+ Template/Probe	5.0	pro Well	5.0	pro Well	5.0	pro Well
= Reaktion gesamt	20.0		25.0		20.0	

2. Proben verteilen (Template je Well separat zugeben)

	Gastro Virus	Gastro Bacteria	Gastro Ecoli	Gastro Virus	Gastro Bacteria	Gastro Ecoli	Gastro Virus	Gastro Bacteria	Gastro Ecoli				
	1	2	3	4	5	6	7	8	9	10	11	12	
A	(-) PTC	(-) PTC	(-) PTC	S240009	S240009	S240008	S240018	S240019	S240016				
B	S240001	S240001	S240001	S240009	S240009	S240009	S240017	S240017	S240017				
C	S240002	S240002	S240002	S240010	S240010	S240010	S240018	S240018	S240018				
D	S240003	S240003	S240003	S240011	S240011	S240011	S240019	S240019	S240019				
E	S240004	S240004	S240004	S240012	S240012	S240012	S240009	S240009	S240009				
F	S240005	S240005	S240005	S240013	S240013	S240013	S240021	S240021	S240021				
G	S240006	S240006	S240006	S240014	S240014	S240014	S240022	S240022	S240022				
H	S240007	S240007	S240007	S240015	S240015	S240015	(-) MTC	(-) MTC	(-) MTC				

260627_Gastro_01

Seite 1 von 1

EXAMPLE OUTPUT

The pipetting scheme as a PDF – MasterMix table and plate overview, colour-matched to the plate setup and to the LightCycler PRO programming. One click, ready to print for the bench.



Download WellDone!

welldone-prompter.github.io/WellDone

Always the latest version – no installation, no administrator rights
Detailed manual (German / English) as a PDF on the same page



SISTER APP

PROMpter

Turns the LightCycler PRO results into a PDF report & Excel automatically after the run.

Sharing is explicitly encouraged – no sale

Feedback: well.done.pcr.setup@gmail.com

Private project, built with Claude Code

Important note: WellDone! is a helper tool, **not a medical device** and not validated. Check the plate setup and the pipetting scheme before use; the professional responsibility for every run stays with the user. No liability is accepted for damage or incorrect results.

"Roche" and "LightCycler" are trademarks of their respective owners. WellDone! is an independent, private project, neither reviewed nor endorsed by Roche. Open-source libraries used: PySide6/Qt, ReportLab.